

01 Vertical SHR With Two-Handle

Fast hair removal, fire up to 10 times/second (10HZ), as fast as Diode Laser

1. Lamp imported from UK, 400,000~500,000 shots per handle can be reached
2. Apply for all hair colors and skin types (I, II, III, IV, V, VI)
3. There are Five filters for your choosing:

- 480nm: vascular therapy
- 480nm: acne therapy
- 530nm: pigmentation therapy
- 530nm: skin rejuvenation
- 590nm: wrinkle removal
- 640nm: hair removal (light colors)
- 690nm: hair removal (deep colors)

TREATMENT HANDLE

Safe and efficient



PRINCIPLE OF TREATMENT

Safe and efficient



Power	3000W
Frequency	1Hz ~ 10Hz
Energy range	1~50 J/cm2
Spot size	15*50mm & 12*30mm
Life time	400,000~500,000 shots per handle

02 Portable SHR With Two-Handle

Fast hair removal, fire up to 10 times/second (10HZ), as fast as Laser Diode.

1. Lamp imported from UK, 400,000~500,000 shots can be reached
2. Apply for all hair colors and skin types (I, II, III, IV, V, VI)
3. There are Five filters for your choosing:

- 480nm: vascular therapy
- 480nm: acne therapy
- 530nm: pigmentation therapy
- 530nm: skin rejuvenation
- 590nm: wrinkle removal
- 640nm: hair removal (light colors)
- 690nm: hair removal (deep colors)

TREATMENT HANDLE

Safe and efficient



PRINCIPLE OF TREATMENT

Treatment process and effect



Power	2500W
Frequency	1Hz ~ 10Hz
Energy range	1~ 50J/cm2
Spot size	15*50mm & 12*30mm
Life time	400,000~500,000 shots
Temperature of handle	-5°C ~ 5°C

03 600W Diode Laser

755nm, 808nm, 940nm, 1064nm, four wavelengths integrated in the same device

1. Two operating systems optional, 808nm system or 755+ 808+940+1064nm system
2. 12 best-quality laser bars imported from "USA-COHERENT"
3. Painless and Permanent result
4. Treatment Handle (with 5 removable heads) or (with adjustable UI screen) optional

TREATMENT HANDLE

Safe and efficient



Power	3000W
Wavelength	755nm, 808nm, 940nm, 1064nm
Laser intensity	5% ~ 100%
Frequency	1Hz ~ 10Hz
Fluence	0~120J/cm2
Pulse duration	5ms ~ 400ms

04 4D Diode Laser

755nm, 808nm, 940nm, 1064nm, four wavelengths integrated in the same device

1. Fast hair removal for all skin types and all hair colors
2. Painless and permanent treatment result
3. Cooling system: semiconductor+ wind+sapphire+TEC
4. Treatment Handle (with 5 removable heads) or (with adjustable UI screen) optional



Power	2500W
Wavelength	755nm, 808nm, 940nm, 1064nm
Laser intensity	5% ~ 100%
Frequency	1Hz ~ 10Hz
Fluence	0~120J/cm2
Pulse duration	5ms ~ 400ms

05 1200W Diode Laser

Permanent hair removal, faster, safer, more effective

1. Two operating systems optional, 808nm system or 755+ 808+940+1064nm system
2. 12 best-quality laser bars imported from "USA-COHERENT"
3. Painless and Permanent result
4. Updated treatment handle with 5 removable tips and adjustable UI screen

TREATMENT HANDLE

Safe and efficient



Output power	1200W
Frequency	1Hz~10Hz (adjustable)
Pulse duration	10ms ~ 200ms
Energy	10 ~ 120J/cm2
Spot size	Φ8, 5*5, 10*10, 12*20, 20*40 (optional)
Screen	10.4" touch LCD + 2" touch LCD

06 755nm + 1064nm Long Pulse Laser

755nm & 1064nm long pulse laser for all skin types, including tanned skin.

1. Solid-state laser and optical Fiber coupling technology, good beam quality
2. Single pulse with high output power, more effective
3. 2~ 20 mm adjustable spot diameter (without changing handle) for treating all body parts.
4. For all skin types including tanned or darker skin

PRINCIPLE OF TREATMENT

Treatment process and effectiveness



BEFORE TREATMENT

Process of treatment



AFTER TREATMENT

Process of treatment



Wavelength	755nm, 1064nm
Output method	Optical Fiber
Output energy	53J (755nm) 80J (1064nm)
Repetition frequency	1-10HZ
Spot diameter	2-20mm (adjustable)
Pulse width	3-200ms

07 Super Picosecond Laser

PicoSure Laser is effective for tattoo, nevus of Ota, coffee spot, freckle chloasma, skin whitening, etc.

1. South Korea imported 7 joint articulated arm
2. 2-10mm adjustable spot size
3. Single pulse(SP), double pulse(PTP), long pulse(LP) optional
4. Solid laser seed source+ MOPA amplification technology
5. Flat-top hat beam mode, Uniform energy density

PRODUCT DETAILS:



Wavelength	1064nm, 532nm (Standard),585nm,650nm (Optional)
Output energy	1300mj(1064nm), 550mj(532nm)
Pulse duration	600PS@1064nm, 500PS@532nm
Spot size	2mm - 10mm
Frequency	1Hz ~ 10Hz

08 Portable Picosecond Laser

PicoSure Laser is effective for tattoo, nevus of Ota, coffee spot, freckle chloasma, skin whitening, etc.

1. South Korea imported 7 joint articulated arm
2. 2-10mm adjustable spot size
3. Single pulse(SP), double pulse(PTP), long pulse(LP) optional
4. Solid laser seed source+ MOPA amplification technology
5. Flat-top hat beam mode, Uniform energy density

PRODUCT DETAILS:



Wavelength	1064nm, 532nm (Standard),585nm,650nm (Optional)
Output energy	1300mj(1064nm), 550mj(532nm)
Pulse duration	600PS@1064nm, 500PS@532nm
Spot size	2mm - 10mm
Frequency	1Hz ~ 10Hz
Voltage	220V/50HZ, 110V/60HZ

09 Diode + SHR + Nd:YAG Laser

3 Handpieces /12 Functions /10 Filters/Safe & Effective

1. Fast Hair removal, Tattoos, Birthmark, Pigmentation, Acne, Vascular, etc
2. UK imported Lamps, effective treatment energy
3. 3 handpieces with 10 filters optional
4. For all skin types including tanned or darker skin
5. Painless and immediate result, no harm, no scar

▶ TREATMENT HANDLE



▶ BEFORE AND AFTER TREATMENT



Power 3000W

Quantity of handle 3

Screen 10.4" color touch LED display

Frequency 1Hz ~10Hz

Cooling system Water+wind+semiconductor+cooling gel

Voltage 220V/50HZ, 110V/60HZ

10 Diode + ND: YAG Laser

2 Handpieces / 9 Functions/ Safe & Effective

1. Fast Hair removal, Tattoos, Birthmark, Pigmentation, Hollywood peel, etc
2. UK imported Lamps, effective treatment energy
3. The lightest handpieces in the world with high power and big spot size, could keep working for 36 hours.
4. For all skin types including tanned or darker skin
5. Painless and immediate result, no harm, no scar

▶ APPLICATIONS

- ✓ Hair Removal
- ✓ Skin Rejuvenation
- ✓ Wrinkle Removal
- ✓ Hollywood Peel
- ✓ Birthmark
- ✓ Onychomycosis
- ✓ Tattoo Removal
- ✓ Pigmentation Therapy
- ✓ Nevus of Ota

▶ PRINCIPLE OF TREATMENT



Power 2500W

Quantity of handle 2

Frequency 1Hz ~ 10Hz

Wavelength 755+808+1064nm(Diode) 532+1064+1320nm(Nd:yag)

Spot size 15mm*15mm or 20*20cm(Diode) 6mm (Nd:yag)

Cooling system Water+Wind+Semiconductor

11 Q-Switch ND:YAG Laser

3 types of laser (532nm, 1064nm, 1320nm) optional, fast tattoo removal

1. 532nm, 1064nm, 1320nm, 3 types of laser.
2. Suitable for all skin types and tattoos.
3. Maximum power can reach 1600mj. Similar products are only 1000mj.
4. This laser can fire up to 10 shots/second (1-10Hz). The treatment is much faster. But the normal laser is only 1-5Hz or 1-6Hz.

APPLICATIONS:



TREATMENT HANDLE



Power	2500W
Pulse duration	8ns
Energy	20MJ ~ 1600MJ
Frequency	1Hz ~ 10Hz (adjustable)
Wavelength	532nm, 1064nm, 1320nm
Voltage	220V/50HZ, 110V/60HZ
Package size	64*53*54cm, 35.5kg

12 4D HIFU + Vaginal HIFU

7 treatment cartridges for multiple treatment needs, such as wrinkle removal, skin rejuvenation, face lift, body slimming, etc.



Maximum power	1000W
4D HIFU cartridges	1.5mm, 3.0mm, 4.5mm, 6.0mm, 13mm
Vaginal HIFU cartridges	3.0mm, 4.5mm
Width	0mm ~ 11mm
Length	5mm ~ 25mm
Frequency:	50-60Hz

13 CO2 Fractional Laser

CO2 Laser for gynecology / dermatology/ plastic surgery with 60W RF laser tube

1. Best brand 'USA-Synrad' RF Laser tube
2. South Korea imported 7 joint articulated arm with good lens
3. Germany imported Mini-Stepping-Motors
4. Fractional mode, Pulse mode, Gynae mode, 3 working modes
5. High quality power supply, ensuring effective power

PRODUCT DETAILS:




Power	1~60W
Scan scope	Max. 20mm * 20mm
Spot diameter	100um ~ 2000um
Treatment mode	Fractional mode, Pulse mode, Gynae mode
Scan graphics	Circle, Square, Triangle, Rectangle, Hexagon, Straight line
Beam delivery	7-joint articulated arm

14 CO2 Fractional Laser

CO2 Laser for gynecology / dermatology/ plastic surgery with 40W RF laser tube

1. Best brand 'USA-Synrad' RF Laser tube
2. South Korea imported 7 joint articulated arm with good lens
3. Germany imported Mini-Stepping-Motors
4. Fractional mode, Pulse mode, Gynae mode, 3 working modes
5. High quality power supply, ensuring effective power

PRODUCT DETAILS:



PRODUCT DETAILS:



Power	1~40W
Scan scope	Max. 20mm * 20mm
Spot diameter	100um ~ 2000um
Mode of treatment	Fractional mode, Pulse mode, Gynae mode
Scan graphics	Circle, Square, Triangle, Rectangle, Hexagon, Straight line
Beam delivery	7-joint articulated arm


15 1927nm Thulium laser

Equipped with fractional handle, ablative and non-ablative double effects


1. Dual ablative and non-ablative effect, safe and short recovery time
2. Low maintenance cost owing to Non-consumable and long life time of Thulium fiber
3. Various sizes of tips and roller tips optional
4. Treat diverse indications by separate setting of "Pulse duration" and "Power"
5. For superficial epidermal pigmented lesions, improving skin texture

PRODUCT DETAILS:

non-ablative
Er.Glass 1550nm



sub-ablative
Thulium 1927nm



ablative
CO2 10600nm



Wavelength	1927nm
Laser power	10W
Beam delivery	Fiber
Spot size	0 ~ 200um
Pulse width	0.1us - 5ms
Scan size	1*1mm ~ 20*20mm

16 980nm 30W Diode Laser


For Liposuction, Vascular, Gynecology, Podiatry, Dermatology, Dentistry, Otorinology

1. Optical Fiber coupling technology, good beam quality.
2. Auto Calibration of the output power at the end of the fibers.
3. One machine with multiple functions.
4. 7 different accessories optional, meet your different treatment needs.
5. 10.4-inch color touch big screen.

PRODUCT DETAILS:



PRINCAIPLE OF TREATMENT



Power	30W
Wavelength	980nm
Output method	Optical Fiber
Operation mode	CW, Single, Repeat
Fibers	400um 600um, SMA905 type connector
Indicator light	650nm red light, power <5mw
Voltage	220V/50HZ, 110V/60HZ

17 4-Handle EMSculptor

THE MORE MUSCLE, THE FASTER FAT BURNS

1. Electromagnetic Wave + Radio Frequency technologies work together for greater energy and efficiency
2. Lying for 30 minutes \approx 5.5 hours of exercise
3. 4 paddles, suitable for treatment of multiple parts of the body
4. Build muscle and burn fat with no surgery and no downtime



Power consumption	2.3KVA
Output intensity	0-14 Tesla
Treatment handle	4 (A, B, C, D)
Contraction (30min)	More than 20000 times
Technology	Electromagnetic Wave + Radio Frequency
Pulse width	300 μ s

18 2-Handle EMSculptor

THE MORE MUSCLE, THE FASTER FAT BURNS

1. Electromagnetic Wave + Radio Frequency technologies work together for greater energy and efficiency
2. Lying for 30 minutes \approx 5.5 hours of exercise
3. 2 paddles, suitable for treatment of multiple parts of the body
4. Build muscle and burn fat with no surgery and no downtime



Power consumption	2.3KVA
Output intensity	0-14 Tesla
Treatment handle	2 (A, B)
Contraction (30min)	More than 20000 times
Technology	Electromagnetic Wave + Radio Frequency
Pulse width	300 μ s

19 2 in 1 Cryo + EMSculptor

Cryo + EMSculptor, lose weight while building muscle

Cryolipolysis is a non-traumatic fat elimination method. Freezing fat leads to fat cell death. At the same time, combined with EMSculptor, it can better shape and increase muscle in a non-invasive way.

PRODUCT DETAILS:



PRODUCT APPLICATIONS



Power consumption 2.3KVA

Output intensity 0-14 Tesla

Treatment handle 4 (A, B, C, D)

Pressure 0.08Mpa


Technology EMSculptor + Cryolipolysis

Voltage 100-220V, 50/60Hz

20 2-Handle Cryolipolysis 360°

Cooling sculpting system for fat reduction and body shaping

1. Two handpieces can work together, shorter treatment time.
2. Two handles can work either independently or simultaneously, more convenient.
3. High-quality refrigeration tablets make sure precise temperature control.
4. 35 ~ 45mins instead of 60mins, save more treatment time.
5. Decrease fat quantity, refuse to rebound.
6. Non-Invasive, Non surgical, No downtime.



PRODUCT DETAILS:



Power 1000W

Temperature -10°C~10°C

Negative pressure 5~80KPa

Handpiece A (big), B (small)

Working mode A, B, A+B

Screen 10.4" color touch LED 3.5" color touch LED*2